

In vitro kinase assay

VD Vladimir Denic HY Houqing YU

Updated date: Sep 1, 2022

 An abbreviated version of this protocol was published in eLIFE in Apr 2022

The peroxisomal exportomer directly inhibits phosphoactivation of the pexophagy receptor Atg36 to suppress pexophagy in yeast

DOI: 10.7554/eLife.74531

Detailed protocol

In vitro kinase assay

Used in Yu *et al.*, *Elife*. 2022 Apr 11;11:e74531.

1. Prepare 20× kinase buffer, aliquot and store at -20 °C.
20× kinase buffer: 1 M Tris-HCl pH 7.5, 2 mM EDTA and 40 mM DTT
2. Prepare fresh 10× ATP regeneration system (ARS) for each assay, keep on ice.
10× ARS: 100 mM MgCl₂, 10 mM ATP, 100 mM creatine phosphate and 2 mg/mL creatine kinase
- 3*. Prepare 4× kinase mixture, keep on ice
0.8 μM Hrr25 (wild type or K38A mutant), 0.8 μCi/μL γ-³²P-ATP, 4× ARS and 4× kinase buffer;
used at 3 μL per 12 μL reaction in PCR tubes
4. Prepare 4× substrate mixture, keep on ice
2 μM Atg36 or 4 μM Atg19, or 2 μM Atg36 and Atg19 diluted with 25 mM HEPES pH 7.5, 150 mM NaCl, and 1 mM DTT
used at 3 μL per 12 μL reaction in PCR tubes
5. Prepare 2× ATPase mixture, keep on ice
400 nM ATPases (wild-type Pex1/6 and Walker B mutant (E832A), wild-type Cdc48 and D2 mutant (E588Q))
used at 6 μL per 12 μL reaction in PCR tubes.
6. Mix 4× substrate with 2× ATPase mixture, pre-incubate the substrate/ATPase mixture on a 37 °C heating block for 5 min
- 7*. Pre-incubate the 4× kinase mixture on a 37 °C heating block for 5 min
- 8*. Add substrate/ATPase into kinase mixture to start the reaction
- 9*. Take 2 μL sample at designed time points (e.g., 0.5, 3, 6, 10 and 15 min), mix the samples immediately with 15 μL 1× SDS-PAGE loading buffer, 95 °C denaturation for 5 min
- 10*. Load the sample (~17 μL) onto SDS-PAGE gels (Novex™ WedgeWell™ 4 to 20%, Tris-Glycine, 1.0 mm, Mini Protein Gels), electrophoresis at 200 V for 55 min
- 11*. Soak the protein gels in gel fixing solution (50% v/v ethanol and 10% v/v acetic acid in water), shake at room temperature for 15 min.
- 12*. Soak the protein gels in gel drying solution (45% v/v ethanol and 5% v/v glycerol in water), shake at room temperature for 15 min.
- 13*. Dry the gels with Model 583 Gel Dryers (Bio-Rad)
- 14*. Expose the dried gels to phosphorimager cassettes overnight
15. Quantify the ³²P-labeled species of Atg36 or Atg19 with FJI software.

*. ³²P related operation. wear appropriate personal protective equipment

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Denic, V. and YU, H. (2022). In vitro kinase assay. Bio-protocol Preprint. [bio-protocol.org/1895](https://doi.org/10.21203/rs.3.rs-1895).
2. Yu, H., Kamber, R. A. and Denic, V. (2022). The peroxisomal exportomer directly inhibits phosphoactivation of the pexophagy receptor Atg36 to suppress pexophagy in yeast. eLIFE. DOI: [10.7554/eLife.74531](https://doi.org/10.7554/eLife.74531)

Copyright: Content may be subjected to copyright.